

ADENOVIRUS CHROMATOGRAPHIC PURIFICATION PROTOCOL

This protocol is adapted from BD Bio Sciences-Clonetech Protocol by the Gene Expression Lab.

This protocol is for use with BD Bio Sciences-Clonetech Systems. For additional technical inquiries, contact Technical Service at 800-662-2566 or www.clonetech.com

ADENOVIRUS AMPLIFICATION
FILTER PURIFICATION
PURIFICATION
VIRUS ELUTION
VIRUS ESTIMATION

BEFORE STARTING THE EXPERIMENT

- Pre-warm all media to RT before use
- Prepare Solutions:
 - 1) Prepare 1X dilution and wash buffer from 5X stock in mill-Q water.
 - 2) Mix the **Benzonase Nuclease** treated filtrate with equal volume of 1X dilution buffer.

Adenovirus Amplification

- 1) Grow 293 cells in 10x150mm plates for 4 days before infection (100% confluent).
 - 2) Infect the 293 cells in 2% media.
 - 3) Incubate the cells at 37°C until Cytopathic effects are complete till cells starts floating.
 - 4) Harvest cells, centrifuge at 3000-RPM, store the cell pellet and supernatant separately.
 - 5) Resuspend the cell pellet in 5 mL of PBS. Lyse cells with 4-5 freeze cycles. Centrifuge and transfer the CVL to a sterile tube.
 - 6) Store the CVL and supernatant in -80C until further use.
-

Filter Purification (Removal of Cell Debris)


Carry out most of the following steps in BL-2 Hood unless indicated otherwise.

- 1) Pool the supernatant and CVL.

- 2) Unwrap a bottle –top filter (from the Kit package) and place in the biosafety hood. Remove the lid from the top of the unit and carefully place the pre-filter disc on top of the 0.45-micron filter.
- 3) Connect the bottle-top filter to a vacuum source.
- 4) With the vacuum off, add 10 mL of 5% media to the pre filter so that it adheres uniformly to the 0.45-micron filter below it (Be careful not to dislodge the pre-filter).
- 5) With vacuum off, carefully pour the pooled supernatant into the top of the filter unit. Being careful to not overfill.
- 6) Turn the vacuum on. When the bottom collection vessel becomes full, disconnect the vacuum line and unscrew the collection vessel.
- 7) Save the filtrate in a sterile flask
- 8) To remove any contaminating cellular DNA, add **Benzonase Nuclease** (Novagen) to the filtrate. The final concentration should be 10 units Benzonase per mL of the filtrate and incubate at 37°C for 30 min.

Purification:

- 1) Collect the following components and place them within easy reach.
 - Tubing assembly and BD adeno X purification filter.
 - 1X wash buffer
 - 1X dilution buffer
 - Sterile PBS
 - 5,30 mL syringes
 - Container for holding biohazardous liquid waste.
- 2) Connect the BD adeno X purification filter and tubing assembly as shown in Figure below.



- 3) Do not place the tubing in the virus-containing solution yet.

- 4) Remove air from the filter and tubing assembly by placing the tube in 50 mL tube containing PBS. Fill a 10 mL syringe with 5-10 mL of sterile PBS.
- 5) Attach the syringe to the filter and push the PBS through the assembly. Place the tubing line in the virus containing solution.
- 6) Pull the virus containing solution through filter at ~20 mL/min. Pull the solution through the filter with the attached 60 mL syringe. When the syringe is full, clamp the tubing assembly, disconnect the syringe, and empty the syringe into a biohazard waste container. Reattach the syringe, unclamp the line, and continue pulling the solution through the filter.
- 7) Repeat this steps until the entire volume has been filtered. Avoid drawing air into the system.
- 8) Wash the adenovirus collected on the membrane by transferring the tube assembly from the virus solution to the 1X wash buffer container.
- 9) Pull the entire volume of the wash buffer through the filter at 20min/mL slowly as described above. Replace the 60 mL syringe with the 3 mL syringe containing elution buffer.

Virus Elution

- 1) Fill a 3 mL syringe with 1X elution buffer. Connect the syringe as described above [to the female port on the filter cartridge,] and remove the tube assembly from the other side.
- 2) Push 1mL of elution buffer through the filter into a sterile 15mL conical centrifuge tube.
- 3) Incubate the filter at room temperature for 5 min , then push the remaining elution buffer through to collect the remaining adenovirus.
- 4) Use air to push any remaining virus from the filter.

Virus Estimation

- 1) Use air to push any remaining virus from the filter.
 - 2) Use elution buffer for blank. Make sure the cuvette is clean outside without any fingerprints.
 - 3) Read O.D at 260nm. Record OD reading. Estimate the amount of virus. (Virus (pfu per mL) = OD at 260 nm x5 x10X10E10 pfu/mL)).
 - 4) Add equal volume of dialysis buffer and store at -70C.
 - 5) Determine the adenoviral titer by BD Adeno-X rapid titer kit.
-